

ELMASONIC MED Verdünnungstabelle



| Gerätetyp | Arbeitsvolumen in ml | Menge in ml | Verdünnung in % | | | | | | |
|--|----------------------|----------------|-----------------|---------|-------|---------|-------|-------|-------|
| | | | 1.0 | 1.5 | 2.0 | 2.5 | 3.0 | 4.0 | 5.0 |
| Elmasonic Med 30 | 1600 | Reiniger in ml | 16 | 24 | 32 | 40 | 48 | 64 | 80 |
| | | Wasser in ml | 1584 | 1576 | 1568 | 1560 | 1552 | 1536 | 1520 |
| Elmasonic Med 60 | 4300 | Reiniger in ml | 43 | 64.5 | 86 | 107.5 | 129 | 172 | 215 |
| | | Wasser in ml | 4257 | 4235.5 | 4214 | 4192.5 | 4171 | 4128 | 4085 |
| Elmasonic Med 100 | 6800 | Reiniger in ml | 68 | 102 | 136 | 170 | 204 | 272 | 340 |
| | | Wasser in ml | 6732 | 6698 | 6664 | 6630 | 6596 | 6528 | 6460 |
| Elmasonic Med 120 | 10000 | Reiniger in ml | 100 | 150 | 200 | 250 | 300 | 400 | 500 |
| | | Wasser in ml | 9900 | 9850 | 9800 | 9750 | 9700 | 9600 | 9500 |
| Elmasonic Med 200 Modularem Korb | 18000 | Reiniger in ml | 180 | 270 | 360 | 450 | 540 | 720 | 900 |
| | | Wasser in ml | 17820 | 17730 | 17640 | 17550 | 17460 | 17280 | 17100 |
| Elmasonic Med 200 mit Sterilisationskorb | 13000 | Reiniger in ml | 130 | 195 | 260 | 325 | 390 | 520 | 650 |
| | | Wasser in ml | 17870 | 17805 | 17740 | 17675 | 17610 | 17480 | 17350 |
| Elmasonic Med 180 | 14200 | Reiniger in ml | 142 | 213 | 284 | 355 | 426 | 568 | 710 |
| | | Wasser in ml | 14058 | 13987 | 13916 | 13845 | 13774 | 13632 | 13490 |
| Elmasonic Med 300 | 22000 | Reiniger in ml | 220 | 330 | 440 | 550 | 660 | 880 | 1100 |
| | | Wasser in ml | 21780 | 21670 | 21560 | 21450 | 21340 | 21120 | 20900 |
| Elmasonic Med 500 | 40300 | Reiniger in ml | 403 | 604.5 | 806 | 1007.5 | 1209 | 1612 | 2015 |
| | | Wasser in ml | 39897 | 39695.5 | 39494 | 39292.5 | 39091 | 38688 | 38285 |
| Elmasonic Med 900 | 63700 | Reiniger in ml | 637 | 955.5 | 1274 | 1592.5 | 1911 | 2548 | 3185 |
| | | Wasser in ml | 63063 | 62744.5 | 62426 | 62107.5 | 61789 | 61152 | 60515 |